


Order form for custom-made

**Compression stockings/panty hose
flat knitted/with seam**



Juzo[®]
Since 1912

Instructions for precise measuring for Juzo® compression stockings/pantyhose

Please note:

- Legs should be free of oedema.
- Measurement should be carried out in the morning since the legs are not or are the least swollen at this time.
- For measuring please use a non-elastic tape measure.
- Measuring should be taken while patient is standing on measuring board (for phlebology and scar treatment).
- All circumferences to be taken with close-fitting tape measure.
- For each circumference you need to take the appropriate length measurement to ensure correct fit as well as the path of compression.
- Length measurements to be taken at inner side of leg/foot.
- Body parts should reach up to middle of waist to be safe from slipping down.

Models:	CCL-AT	CCL-Leg segments	CCL-Body part
	1	1	1
	2	2	2
	3	3	2
	4	4	2

The models are standard combinations – every other combination is possible.

Detailed instructions for precise measuring can be found in the Juzo® flyer “Precise measuring for compression stockings and compression pantyhose”.

For lymphological conditions, please note the following additional instructions:

- Take the length measurements ℓD , ℓF and ℓG (ℓG for pantyhose = up to crotch/ ℓG for stockings = up to natal cleft) as well as the trunk section circumferences while the patient is standing.
- Take all other leg measurements while the patient is lying down.
- Do not measure within skin folds.
- When taking circumference measurements under tension, please do not do this jerkily but always in a manner that is comfortable for the patient.
- We recommend angling the openings of the stockings in order to ensure an optimal pressure distribution and anatomical fit.
- Photographic documentation can be of assistance in difficult anatomical situations.

If you want to order combined supplies please use two separate order forms.

For the selection of the appropriate garment, we recommend the brochures Juzo® Product range and Juzo® Expert Prescription brochure.

Model:

- ☐ AD below knee stockings
☐ AG thigh high stockings
☐ AG thigh high with hip attachment
☐ AT tights
☐ capri pants
☐ bermuda pants

_____ piece _____ pair

Colours¹⁾: ☐ Beige ☐ Platinum ☐ Red
☐ Pink ☐ Blue ☐ Grey ☐ Mocca
☐ Dark blue ☐ Black ☐ Green²⁾
 (compression class 4 only available in colour Beige)

¹⁾ If nothing is noted, we deliver colour Beige

²⁾ Delivery time approx. 10 days

Options (no charge):

Stockings:

- ☒ full knit in the thigh part
 (standard for Expert and Expert Strong)

Tights:

- ☐ box type (standard: diagonal version)
☐ standard border without waist band at "T"
 Gusset: ☐ normal ☐ small
 Mesh gusset: ☐ normal ☐ small
☐ one legged pantyhose

Charge options:

Stockings:

- ☐ overheight
☐ overheight max
☐ adhesive border
☐ lace border
☐ decorative adhesive border
 Hip attachment (please mention "cT")
☐ left ☐ right ☐ "to be worn as a pair"

Tights:

- ☐ with adjustable waist band
☐ at "T" band of the hip attachment
☐ scrotum included
☐ fly
☐ open crotch
☐ body part in the middle of the front with zipper (not with adjustable waist band)
☐ Pressure pads designed by Dr. Rotter
☐ left ☐ interior ☐ exterior
☐ right ☐ interior ☐ exterior

Toe options:

- ☐ open toe (no charge)
☐ closed toe (charge option)
☐ additional instep stitching (charge option)
☐ ball stub/diagonal foot (charge option)

"lA" left _____ cm interior

_____ cm exterior

"lA" right _____ cm interior

_____ cm exterior

From:

Date: _____

Your order no./Commission: _____

Expert
Expert Strong

Compression classes

	1	2	3	4
Expert	<input type="checkbox"/> 3021	<input type="checkbox"/> 3022	<input type="checkbox"/> 3023	<input type="checkbox"/> 3024
Expert Strong	<input type="checkbox"/> 3051	<input type="checkbox"/> 3052	<input type="checkbox"/> 3053	<input type="checkbox"/> 3054

Patient ☐ male
☐ female

Former order no./dd.: _____

Circumferences in cm		Length of body part K-T in cm	
left	right		
cT.....		front	
cH.....		back	
cK.....		Length measurements in cm	
		left	right
cG.....	cG.....	lT	
cF.....	cF.....	lH	
cE.....	cE.....	lG/lK	
cD.....	cD.....	(for panty hose crotch measurements)	
cC.....	cC.....	lF	
cB ¹	cB ¹	lE	
cB.....	cB.....	lD	
cY.....	cY.....	lC	
cA.....	cA.....	lB ¹	
		lB	
		lA Foot length (Back of heel to cA)	
		lZ Total foot length (indication only if desired with closed toe)	

Special requests: _____

(Further options see Juzo® Expert Prescription brochure)

Juzo UK Ltd
South Court
Sharston Rd
Sharston
Manchester M22 4SN
Great Britain

Telephone 01382 826620
Telefax 01382 826641
E-mail: sales@juzo.com

www.juzo.com